

PREPARATION FOR THE CONJUNCTIVAL SWAB

INDICATIONS

The examination is performed in the **Seriato offices** from Monday to Friday during the opening hours of the sampling room and on Saturdays by appointment; **Alzano Lombardo, Piario, Lovere** from Monday to Friday during the opening hours of the sampling room.

Suspend the instillation of ophthalmological drops in the 24 hours prior to sampling. Suspend any antibiotic therapy 4 days before the examination.

COLLECTION MATERIAL

- AMIES Pad with Transport Media (PINK CAP)

WITHDRAWAL METHODS

Collect secretions or samples of tear fluid or mucus from the inferior conjunctiva (surface of the eye) by gentle rotation of the swab.

